


Rapid COVID-19 Antigen Self-Test

Benefits

- ▶ Allow the individuals to collect specimens by themselves or with the assistance of others.
- ▶ Can detect all suspected case specimens or general screening specimens.
- ▶ Simple test procedure using a nasal swab from the front of the nasal cavity instead of the nasopharynx.
- ▶ With single use components, safe & easy operation and reduces the risk of pollution.
- ▶ Receive highly accurate results within 15 minutes.

Kit Contents

- 1/2/3/5/7/10/15/20 Test Device(s)
- 1/2/3/5/7/10/15/20 Sterile Swab(s)
- 1/2/3/5/7/10/15/20 Extraction Tube(s) with Buffer and Tip(s)
- 1-2 Workstation(s) (If applicable)
- 1 Instructions For Use
- 1 Procedure Card




Test Procedure




Step 1.
Wash or sanitize your hands before performing the test.



Step 2.
Blow your nose. Keep 1.5 meters away from bystanders when blowing your nose.



Step 3.
Insert the tube into the workstation. Remove the lid from the top of the tube.




Step 4.
Open the swab package where indicated. Pull the swab out by grasping the plastic end. Do not touch the absorbent swab tip.



Step 5.
Carefully insert the absorbent tip of the swab into your left nostril. Ensure that the entire swab tip is inside your nostril (2–4 cm deep). Do not insert the swab further after you feel resistance.



Step 6.
Roll the swab at least 5 times against the insides of your nostril. Ensure good contact between the swab and the insides of your nostril.



Step 7.
Remove the swab and insert into your right nostril. Repeat Steps 3 and 4.



Step 8.
Remove the swab from your nostril and insert the swab into the prepared tube in the workstation.




Step 9.
Mix well by rolling the swab at least 6 times while pressing the head of the swab against the bottom and sides of the tube.



Step 10.
Start timer. Leave the swab in the tube for 1 minute.



Step 11.
Squeeze the tube several times from the outside. Try to dissolve as much solution from the swab as possible.



Step 12.
Remove the swab and discard in a plastic bag.



Step 13.
Push the tip provided in the kit into the tube and ensure it fits tightly.



Step 14.
Remove the test device from the pouch and lay it on a flat clean surface.




Step 15.
Add 4 drops of the solution into the sample well of the test device by gently squeezing the tube.



Step 16.
Start timer and read results at 15 minutes. It is important to read the results at 15 minutes.


Step 17.
Read your results. There are three possible types of results.

Positive Test Result




If the test device looks like either of the positive result windows as shown on the left, you have a current Covid-19 infection. Please call your doctor or your local health department and make sure you adhere to local guidelines for self-isolation. Re-testing with other test methods such as a PCR test may be required.

Negative Test Result



If the test device looks like the negative result window as shown on the left, no Covid-19 infection could be detected. In a suspected case, repeat the test after 1-2 days since the virus cannot be accurately detected in all phases of an infection. Despite a negative test result, you still have to comply with all applicable rules regarding contact with others and protective measures.

Invalid Test Result



If your test result looks different, meaning there is no line visible or only one line at T, the result is invalid. This may be a result of the test execution, and the test should be repeated. If invalid test results continue, please contact your doctor or a COVID test center.

Ordering Information

Product Description	Specimen	Catalog No.	Format	Kit Size
Rapid COVID-19 Antigen Self-Test	Nasal Swab	GCCOV-502a-H1 ✓	Cassette	1 Test/Kit
Rapid COVID-19 Antigen Self-Test	Nasal Swab	GCCOV-502a-H2 ✓	Cassette	2 Tests/Kit
Rapid COVID-19 Antigen Self-Test	Nasal Swab	GCCOV-502a-H3 ✓	Cassette	3 Tests/Kit
Rapid COVID-19 Antigen Self-Test	Nasal Swab	GCCOV-502a-H5 ✓	Cassette	5 Tests/Kit
Rapid COVID-19 Antigen Self-Test	Nasal Swab	GCCOV-502a-H7 ✓	Cassette	7 Tests/Kit
Rapid COVID-19 Antigen Self-Test	Nasal Swab	GCCOV-502a-H10 ✓	Cassette	10 Tests/Kit
Rapid COVID-19 Antigen Self-Test	Nasal Swab	GCCOV-502a-H15 ✓	Cassette	15 Tests/Kit
Rapid COVID-19 Antigen Self-Test	Nasal Swab	GCCOV-502a-H20 ✓	Cassette	20 Tests/Kit

✓ CE Marked

